

# Product Details

## Product Information

Galilean Loupes have the advantage of being easy to use, with a high depth of field that keeps things in focus while moving. The surgical medical binocular loupes adopt lightweight aluminum strips, which can independently adjust the pupil, make the installed magnifier stable, ideal focus calibration, avoid headache and eye fatigue.

Surgical medical binocular loupes are specially designed for comfort and high performance. The distance during treatment can make the patient feel better and easier. Lister magnifiers can perform according to your requirements. The Galilean Loupes have multiple magnification levels so you can see more spot-like problems and even provide precise treatment.

## Features of surgical medical binocular loupes:

- Suitable for general dentistry, routine surgical procedures, tooth preparations and dental hygiene.
- Asynchronous adjust type
- Multi-layer anti-reflection coating

## Specifications of GL3

Magnification	Working distance	Field of view	Depth of field	Weight
3.5x	340-500mm	60-115mm	65-175mm	65.0g

## Frame Options:

- Sport frame: E sport frame, D sport frame
- TR frame: AM56276, AM54118
- Aluminum frame:AM161, AM503, AM549
- Titanium Aluminum frame:AM238
- Alu-mag frame: AM559, AM630,AM640

## Case Options



Fabric case