

Traus SIP20

Premium Implant Engine

EZ OPEN

80 Ncm

OPTIC

- Maintenance friendly
- Intuitive user interface
- Real time torque & RPM display
- Motor auto-calibration function
- Efficient gear system with the latest technology



ORDER CODE

· SET	SIP20	
· CONTROL BOX	BIP20	· MOTOR MBP10SL
· FOOT CONTROLLER	FUS10	· ANGLE CRB46LN

Micro Motor & Handpiece



MBP10SL

- Optic : Optic
- Weight : 139 g
- Dimension : 127 × ϕ 25.7 mm
- Led Power : 25,000 Lux



CRB46LN

- Optic : Optic
- Max. RPM : 2,000
- Gear Ratio : 20 : 1
- Min. RPM : 15
- Max. Torque : 80 Ncm
- Weight : 85 g

CONTROL BOX

- Input : AC 100-120 V, 220-240 V (50/60Hz)
- Output : 48 VA
- Fuse : 5 A / 2 A
- Maximum Pumps : Max. 90 $\pm$ 20%ml / min
- Dimension : W 285.5 × D 261 × H 121.7 mm
- Weight : 3.42 kg



Maintenance Friendly

Easy to disassemble by hand without any tools for internal cleaning and maintenance

