


Traus B110



- Knee ↔ Foot control select
- Compact design
- Auto-Cruise function
- Overload protection system

SPECIFICATIONS

· Speed	: Max. 50,000 RPM
· Input	: AC 100-120 V, 220-240 V (50/60Hz)
· Output	: 230 W
· Weight	: 2.4 kg



H180

- Slim & Compact type
- Powerful torque
- Less noise and vibration
- Ergonomic design

SPECIFICATIONS

· Speed	: Max. 50,000 RPM
· Torque	: Max. 8.2 N.cm
· Weight	: 202 g (without cord)



ORDER CODE

· SET	TRAUS B110 + H180 + FOOT CONTROLLER
· CONTROL BOX	TRAUS B110
· HANDPIECE	H180
· FOOT CONTROLLER	SPEED