

## Relaxing Music For a Better Emotional Healing

- Water adjust knob.
- Water/Air Switch.



## Multi-mode Music Playback

- Music playing through USB, Bluetooth, local disk storage and other ways can relieve the tension of patients, eliminate the fear of patients and improve patient's satisfaction.



- The sidebox is connected in the middle of the frame.
- Hanged type aluminium light arm.
- Multi-functional tray, USB Charging port integrated on 3-in-1 cup dispenser.



## 90° Rotary Water Sidebox Element

- Detachable sidebox cover.
- The water can pour into the bottles from the top of the sidebox.
- Muffling device for strong suction.
- Water/Air/Electric separate design.



## Dentist Element

- First-aid chair position.
- One-touch self-stop function.
- Clean up chair position.



- Built-in X-ray film viewer.
- Full size acrylic cover.
- 6 -rack design.



## First-class Comfort Experience

- Seamless high quality PU upholstery in Ergonomics design.



- Double-joint headrest.
- Alarm button.
- Ultra-thin wing-shaped metal backrest.



## Cingol 3.0 Smart Operation System

- Seat Decline protection.
- Skidproof handrest.
- Snap design for front and rear cover.
- Emergency alarm.



## Assistant Panel

- Multi-functional control with 4-rack design.
- 180° angle limit to avoid touching the chair backrest.
- Anti-perforation infrared touch panel.



## Operating Lamp

- 3-joint design reflecting type operating lamp.
- USA CREE LEDs, high lumen & CRI.
- Eco-friendly and durable.
- Dual Color ( Yellow & White ) ,Close to natural light.



## 180° Rotary Ceramic Spittoon

- Detachable ceramic spittoon.
- Easy for disinfection and cleaning.
- 180° rotatable spittoon.
- Deepening the design to avoid splashing.




## Multi-function Foot Pedal

- Cup filling, chip blow, LED lamp, chair movement.
- Air/water supply separately.
- Interlock machine chair.



Input voltage: 230V  $\pm$  10% ,50Hz  
Water pressure: 0.2-0.4Mpa/2.0-4.0 bar  
Air pressure: 0.55-0.8Mpa/5.5-8.0 bar  
G.Weight: 255Kg



## Color Options

