

MMPP

MILLIMETER MARKED PAPER POINTS



- MMPP is marked to provide preliminary guidance in measuring the root canal depth before proceeding clinically with the application of the Gutta Percha Points.
- These points are highly absorbent and color coded to identify ISO tip sizes.

Available Sizes

Size	SP-6 0
XX F	MP202-S601
XF	MP202-S602
F	MP202-S603
M	MP202-S604
C	MP202-S605
XC	MP202-S606
XXF XC	MP202-S691

(Pkg. of 200)

Size	SP-6 0
#8	MP201-S60 1
#10	MP201-S60 2
#15	MP201-S60 3
#20	MP201-S60 4
#25	MP201-S60 5
#30	MP201-S60 6
#35	MP201-S60 7
#40	MP201-S60 8
#45	MP201-S60 9
#50	MP201-S61 0
#55	MP201-S61 1
#60	MP201-S61 2
#70	MP201-S61 3
#80	MP201-S61 4
#15/4 0	MP201-S691
#45/8 0	MP201-S692
#90/14 0	MP201-S693

Package Information

- P-33 Cell Package
Points remain sterile until opened
Box dimension: 78×19×50mm(W×H×D)
- SP-60 Slide Package
One hand slide open lid
Box dimension: 70×13×45mm(W×H×D)

Size vs Color Chart

8	10		
15	45	90	XXF
20	50	100	XF
25	55	110	F
30	60	120	M
35	70	130	C
40	80	140	XC

Contents

Size	P-33	SP-60
#08-80	200p s	200p s
#90-140	200p s	120p s

COLOR CODED PAPER POINTS

- Made from specially formulated highly absorbent papers.
- Color Coded.
- Hand Jig rolled to ensure Superior absorbency.



Package Information

- P-60 Spill Proof Package
6 compartments plastic package
Box dimension:
71×13×50mm(W×H×D)

Available Sizes

Size	P-60
XX F	202-60 1
XF	202-60 2
F	202-60 3
M	202-60 4
C	202-60 5
XC	202-60 6
XXF XC	202-69 1

(Pkg. of 200)

Size	P-60
#8	201-60 1
#10	201-60 2
#15	201-60 3
#20	201-60 4
#25	201-60 5
#30	201-60 6
#35	201-60 7
#40	201-60 8
#45	201-60 9
#50	201-61 0
#55	201-61 1
#60	201-61 2
#70	201-61 3
#80	201-61 4
#15/4 0	201-69 1
#45/8 0	201-69 2
#90/14 0	201-69 3

Size vs Color Chart

8	10		
15	45	90	XXF
20	50	100	XF
25	55	110	F
30	60	120	M
35	70	130	C
40	80	140	XC

Contents

Size	P-60
#08-8 0	200p s
#90-14 0	120p s

DIA-PROISO .04/.06 PLUS

MILLIMETER MARKED SPECIAL TAPER PAPER POINTS

- Millimeter marked Special Taper Absorbency Paper Points.
- Designed to be used for drying and cleaning the canals prepared with **Profile**®, **Hero**® & **K3**® files.
- Available in .04 & .06 taper.



Available Sizes

Size	Dia-Pro ISO .04 Plus	Dia-Pro ISO .06 Plus
	SP-60	SP-60
#15	MP258-S603	MP259-S603
#20	MP258-S604	MP259-S604
#25	MP258-S605	MP259-S605
#30	MP258-S606	MP259-S606
#35	MP258-S607	MP259-S607
#40	MP258-S608	MP259-S608
#45	MP258-S609	MP259-S609
#50	MP258-S610	MP259-S610
#55	MP258-S611	MP259-S611
#60	MP258-S612	MP259-S612
#70	MP258-S613	MP259-S613
#80	MP258-S614	MP259-S614
#15/40	MP258-S691	MP259-S691
#45/80	MP258-S692	MP259-S692

(Pkg. of 100)

DIA-PROISO .04/.06 PLUS IN THE CELL PACKAGE

DIA-PROISO .04/.06 PLUS IN THE CELL PACKAGE



- Points remain sterile until opened

Cell pkg	Dia-Pro ISO .04 Plus	Dia-Pro ISO .06 Plus
#15	MP 258-703	MP 259-703
#20	MP 258-704	MP 259-704
#25	MP 258-705	MP 259-705
#30	MP 258-706	MP 259-706
#35	MP 258-707	MP 259-707
#40	MP 258-708	MP 259-708
#15/40	MP 258-791	MP 259-791

(Pkg. of 120)

DIA-PROT

DIA-PROT



- Designed to be used to complement **DIA-PT**™ File, **Dia-X** File and **Protaper**® Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#F1	MP250-S601
#F2	MP250-S602
#F3	MP250-S603
#F4	MP250-S604
#F5	MP250-S605
#F1/F3	MP 250-S691
#F4/F5	MP 250-S692
#F1/F5	MP 250-S693

(Pkg. of 100)

DIA-PROT NEXT

DIA-PROT NEXT




- Designed to be used to complement **Protaper**® Next Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#X2	MP267-S602
#X3	MP267-S603
#X4	MP267-S604
#X5	MP267-S605
#X2/X5	MP267-S691
#X2/X3	MP267-S692

(Pkg. of 100)

DIA-PRO W



- Designed to be used to complement **Wave One**® Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	MP265-601
#Primary	MP265-602
#Large	MP265-603
#Asst(Small/Primary/Large)	MP265-691

(Pkg. of 100)

DIA-PRO W GOLD



- Designed to be used to complement **Wave One**® Gold Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	MP268-S601
#Primary	MP268-S602
#Medium	MP268-S603
#Large	MP268-S604
#Asst(Small/Primary/Medium/Large)	MP268-S691

(Pkg. of 100)

DIA-PRO R



- Designed to be used to complement **RECIPROC**® Files.
- Millimeter Marked Points

Size	ITEM No.
#25	MP264-605
#40	MP264-608
#50	MP264-610
#Asst(25/40/50)	MP264-691

(Pkg. of 100)