



STERI. IRRIGATION TIPS-SIDE OPEN TYPE

- Disposable Endodontic & Periodontal Irrigation Tips with Safety.
- Unique design with a side-port for superior effectiveness and patient comfort in irrigating root canals and periodontal pockets.
- Pre-sterilized by E.O Gas.
- Use with Luer Lock Syringes.
- Side-Port and Closed-end Irrigation tips.



PERIODONTAL POCKETS

- Side-port design make an effective irrigate pockets and closed, rounded tip can reduce the risk of tissue damage.



ROOT CANALS

- Side-port dispersal can create an upward flushing motion for the complete root canal preparations and cleanse the gingival sulci fully without damaging tissue.

size	Color	ITEM No.
31G (0.25 x 25mm)	Clear White	951-101A
30G (0.30 x 25mm)	Yellow	951-102A
28G (0.36 x 25mm)	Blue Green	951-103A

(Pkg. of 50)

