

LAFOMED



8L

12L



18L

23L

| Program         | 121℃<br>solid | 121℃<br>Universal | 134℃<br>solid | 134℃<br>Universal | BD/Helixa<br>Test | Vacuum<br>Test |
|-----------------|---------------|-------------------|---------------|-------------------|-------------------|----------------|
| Vacuum Times    | 1             | 3                 | 1             | 3                 | 3                 | 3              |
| Temperature (℃) | 121           | 121               | 134           | 134               | 134               | /              |
| Pressure (KPA)  | 110           | 110               | 210           | 210               | 210               | -80            |
| Time (Min.)     | 20            | 20                | 4             | 4                 | 3.5               | 15             |

| Standard 8L  | LF-8L-E          | Standard 12L | LF-12L-E         | Standard 18L | LF-18L-E         | Standard 23L | LF-23L-E         |
|--------------|------------------|--------------|------------------|--------------|------------------|--------------|------------------|
| Capacity     | 8L               | Capacity     | 12L              | Capacity     | 18L              | Capacity     | 23L              |
| Chamber Size | Φ17*32cm         | Chamber Size | Φ20*36cm         | Chamber Size | Φ25*35cm         | Chamber Size | Φ25*45cm         |
| Display      | LED              | Display      | LED              | Display      | LED              | Display      | LED              |
| Lock         | Double lock      | Lock         | Double lock      | Lock         | Double lock      | Lock         | Double lock      |
| Product Size | 63*43*36cm       | Product Size | 63*43*36cm       | Product Size | 62*50*43cm       | Product Size | 66*50*43cm       |
| G.W./N.W.    | 52/45KGS         | G.W./N.W.    | 55/47KGS         | G.W./N.W.    | 61/53KGS         | G.W./N.W.    | 63/55KGS         |
| Power        | 220V/50Hz 1.95KW | Power        | 220V/50Hz 1.95KW | Power        | 220V/50Hz 1.95KW | Power        | 220V/50Hz 1.95KW |
| Output       | Printer          | Output       | Printer          | Output       | Printer          | Output       | Printer          |

LAFOMED

Ningbo Haishu Life Medical Technology Co.,ltd.

No.1 Jinghui Rd., Hengjie Town, Haishu, Ningbo, Zhejiang.  
315181, China PRC

Tel: 0086 574 8700 6890 Fax: 0086 574 8828 3099

E-mail: sales@lafomed.com



LAFOMED



AUTOCLAVE  
STEAM STERILIZER

STANDARD LINE

Prevent Infection. Create Safety.

Standard Line



**B Type**

In compliance with international standards



Built-in printer



**Built-in Printer**

Easy & safe data saving



**Physical Buttons**

Easy operation & longer service life



**4 Preset Program, 2 Test Program**

Suitable for all kinds of medical instruments



**Reliable Locking System**

- Double lock (mechanical and electrical)
- Two position mechanical lock design