


handy™



HDI-220C

Auto-focus(fixed) Intraoral Camera

CE

**Clear Images**  
High-definition image quality with less than 5% distortion helps to show the details perfectly.

**Mechanical Button**  
Mechanical buttons with easy-to-use features.

**Metal Housing**  
Oxide aluminum alloy shell, easy to clean and more durable.

**Professional Lens**  
Professional dental camera lens with longer service life and stronger anti-aging ability.

**Natural-Color LED**  
Six natural-color LED lights, restore the true color of the intraoral images.

**Image Sensor**  
1/3" HD CMOS

**Signal Output**  
USB 2.0

**Resolution**  
1280 x 720 (720p)

**Frame Rate**  
30fps@720p

**Focal Length**  
5mm~35mm

**Lighting**  
6 LEDs

**USB Length**  
2m

**Angle of View**  
≥80°

**Operating System**  
Win7/Win10(32bit&64bit)

**Handy Medical Equipment CO., Ltd**



Tel: +86 021-3616 1602  
Website: [www.handydental.cn](http://www.handydental.cn)  
E-mail: [info@handycreate.com](mailto:info@handycreate.com)

**Scale 1:1**  
**USB Cable Connection**

