

BEAUTIFIL-Bulk Flowable

CLICK / SCAN



FOR MORE INFO

BIOACTIVE BULK FLOWABLE COMPOSITE


Designed for convenient dentin replacement in posterior restorations with self-levelling flowability ...

- **Seamless Adaptation** to cavity walls with curing depth upto 4mm
- **Simplified Inventory** with just 2 shades
- **Bioactive Benefits** to prevent secondary caries

Formulated with a unique balance of filler- resin chemistry, which includes new monomers and a filler complex with the patented S-PRG filler technology to provide, deep cure with minimal polymerization shrinkage stress and optimal physical properties for ensure lasting restorations.



Available in 2 Shades – 2.4gm syringe



CLINICAL APPLICATIONS

- Strong base or liner for Class I and II cavities
- Small posterior restorations
- Undercut block out
- Minor repair of defects on direct resin restorations
- Restoration of deciduous teeth

FEATURES

- ✓ **Ideal flowability and seamless adaptation**
 - Up to 4mm deep cure in a single application
 - Fast & efficient dentin replacement
 - Eliminates voids
- ✓ **Unique Syringe Design**
 - 360° multidirectional grip
 - Plunger designed to minimize wastage with pull-back mechanism
- ✓ **Mimics the shade of natural posterior teeth**
 - 2 versatile shades that blends well
 - Stable shades pre & post cure
- ✓ **Bioactive to protect against secondary caries**
 - Minimize plaque formation, reduce cariogenic bacteria & neutralize sugar acids
 - Fluoride release & recharge
- ✓ **Long lasting restorations**
 - Strong and durable
 - Minimize polymerization shrinkage
 - Increased wear resistance

PROPERTIES

- Flexural Strength: 117 Mpa
- Filler load: 73 wt%
- Radiopaque

RELATED PRODUCTS

- Beautifil Bulk, Beautifil II, Beautifil II LS
BeutiBond / FL Bond II



Official Partner

